



NOT ALL BICARB SYSTEMS ARE CREATED EQUAL CHECK OUT THE FEATURES BEFORE YOU BUY

Bicarb Mix & Distribution System Features	CMS Requirement	AmeriWater
<i>Separate mix & distribution tanks</i>		<i>Option</i>
<i>Cone bottom tanks w/drain at bottom of cone</i>	Yes	Yes
<i>High level & low level alarms</i>	Yes	Yes
<i>Translucent tanks – no site tubes</i>	Yes	Yes
<i>Tight fitting lid on mix & distribution tanks. No opening for paddle type mixer</i>	Yes	Yes
<i>Mix & distribution tank air vent, 0.2um hydrophobic air vent filter</i>	Yes	Yes
<i>Spray ball in mix tank to prevent microbiological growth</i>	Yes	Yes
<i>Mix tank air vent, 0.2um hydrophobic air vent filter</i>	Yes	Yes
<i>Programmable auto fill timer on mix tank</i>		Yes
<i>Manual control for fill on mix tank</i>		Yes
<i>Separate mix & distribution pumps to allow mix of new batch while distributing current batch.</i>		Yes
<i>Simple method to empty mix/storage tank at end of day and rinse system & loop</i>	Recommended	Yes
<i>Method to avoid over mixing bicarbonate (mixer timer)</i>	Yes	Yes

STANDARDS:

- FDA 510K / ISO 13485 Registered Medical Device
- Health Canada Medical Device License 81988



Automatic Dual Tank



Single Tank Mix & Distribution



Manual Mix & Dispense